

NIVA

Lumbar PEEK Spacer System
Posterior Lumbar Interbody Fusion



ASAM

Health is wealth, leave it to us



Features

- ▶ Bulleted nose for easier implantation especially in very degenerated discs
 - ▶ Clamping mechanism for safe and easy connection with the inserter
 - ▶ Wide range of sizes for an individual patient treatment as for example 1 mm increments in height
 - ▶ Optimized ratio between contact area and opening

Introduction

The material used is biocompatible PEEK, PEEK stands for PolyEtherEtherKetone. PEEK polymer complies with ISO 10993-1, USP Class VI and ASTM F2026 for use as a medical implant material.

The use of PEEK as an orthopedic device material enjoys increased popularity in recent years due to the material's unique combination of characteristics. It's properties include radiolucency, high mechanical strength, biocompatibility and compatibility with standard sterilization methods.


The intrinsic radioscopic transparency of the material provides permeability on X-rays and CT scans, allowing to visualize bone growth adjacent to the implant. This allows quick and simple assessment of the bone structure and progress towards bone fusion.





ASAM

Health is wealth. leave it to us

NIVA



NIVA PEEK Cage-ordering information



Ref No.	Height	Width	Length	Angle	Ref No.	Height	Width	Length	Angle
15-092207	7	9	22	0°	16-092207	7	9	22	5°
15-092208	8	9	22	0°	16-092208	8	9	22	5°
15-092209	9	9	22	0°	16-092209	9	9	22	5°
15-092210	10	9	22	0°	16-092210	10	9	22	5°
15-092211	11	9	22	0°	16-092211	11	9	22	5°
15-092212	12	9	22	0°	16-092212	12	9	22	5°
15-092213	13	9	22	0°	16-092213	13	9	22	5°
15-112507	7	11	25	0°	16-112507	7	11	25	5°
15-112508	8	11	25	0°	16-112508	8	11	25	5°
15-112509	9	11	25	0°	16-112509	9	11	25	5°
15-112510	10	11	25	0°	16-112510	10	11	25	5°
15-112511	11	11	25	0°	16-112511	11	11	25	5°
15-112512	12	11	25	0°	16-112512	12	11	25	5°
15-112513	13	11	25	0°	16-112513	13	11	25	5°

Please visit [www. asammedical.com](http://www.asammedical.com) for full information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization.



 No. 104, Between 5 & 7th Talash
St, Kashmar Industrial Estate, Kashmar,
Khorasan Razavi Province, IRAN

www.asammedical.com
info@asammedical.com
Document No: ASM-PC-QA-003
Date: 2023.01.07