



RENUS

Lumbar Posterior
Fixation System

Introduction

The RENUS Spinal Fixation System is a comprehensive implant designed for degenerative and deformity correction procedures in posterior thoracolumbar fusion applications. The system includes Mono & Polyaxial pedicle screw, rod, cross-link and hooks to single level and multi-level to provide stabilization fusion. The system can be rigidly locked into a variety of configurations, with each construct being customized to the patient's anatomy.




Features

- Color code by screw diameter
- Wide range of screw length
- Reduction screws provide controlled and gradual anatomical reduction and simply break them off to remove
- Mono & polyaxial pedicle screw
- Thread design allows for increasing bone stability
- The design allows up to 40° angulation, offering easy connection and greater flexibility
- The cutting flute improves bone insertion and self-tapping
- The nut screw thread design for securing between rod and screw

Reduction screw provide controlled and gradual anatomical reduction and simply break off to remove it

Cutting flute improves bone insertion



Color Coded

Up to 50° angulation offering easy connection and greater flexibility

Thread design increasing bone fixation stability



Nut screw thread design for securing connection between rod and screw

RENUS-ordering information



Polyaxial Pedicle Screw, Long Arm

Ref No.	Diameter × Length
30-135550	50 × 5.5 mm
30-136530	30 × 6.5 mm
30-136535	35 × 6.5 mm
30-136540	40 × 6.5 mm
30-136545	45 × 6.5 mm
30-136550	50 × 6.5 mm
30-137540	40 × 7.5 mm
30-137545	45 × 7.5 mm
30-137550	50 × 7.5 mm
30-137555	55 × 7.5 mm

Ref No.	Diameter × Length
30-134525	25 × 4.5 mm
30-134530	30 × 4.5 mm
30-134535	35 × 4.5 mm
30-134540	40 × 4.5 mm
30-134545	45 × 4.5 mm
30-135525	25 × 5.5 mm
30-135530	30 × 5.5 mm
30-135535	35 × 5.5 mm
30-135540	40 × 5.5 mm
30-135545	45 × 5.5 mm

Nut Screw Lumbar



Ref No.	Size
30-141000	10 mm

Lumbar Cross Link



Ref No.	Size
30-185000	50 mm
30-186000	60 mm
30-187000	70 mm
30-188000	80 mm

Lumbar Straight Rod



Ref No.	Size
30-155505	50 × 5.5 mm
30-155506	60 × 5.5 mm
30-155507	70 × 5.5 mm
30-155508	80 × 5.5 mm
30-155509	90 × 5.5 mm
30-155510	100 × 5.5 mm
30-155511	110 × 5.5 mm
30-155512	120 × 5.5 mm
30-155514	140 × 5.5 mm
30-155516	160 × 5.5 mm
30-155518	180 × 5.5 mm
30-155520	200 × 5.5 mm
30-155530	300 × 5.5 mm
30-155540	400 × 5.5 mm
30-155550	500 × 5.5 mm


Please visit asammedical.com for full information on indications for use, contraindications, warnings, precautions, adverse reactions and sterilization



ASAM

Health is wealth. leave it to us



 No. 104, Between 5 & 7th
Talash St, Kashmar Industrial
Estate, Kashmar, Khorasan Razavi
Province, IRAN

www.asammedical.com
info@asammedical.com
Document No: ASM-PC-QA-004
Date: 2023.04.13